

STtuve® III



IF THE **COLOUR CHANGES** **BACTERIA ARE PRESENT!**

Quick Results in 12 to 24 hours

Accurate Selective liquid media ensures optimal bacterial growth

Simple Take a swab, incubate and observe the colour change

STtuve III is so easy and convenient, anyone can use it

- **Implementation and maintenance of HACCP**
- **Minimise the risk of food poisoning**
 - **Single use, self contained product**
 - **18 month shelf life**



If you are concerned about:

- **The presence of potentially harmful bacteria?**
- **Possible bacterial contamination of your products?**
- **Your ability to respond quickly to bacterial contamination?**
- **Having to wait days for results?**

Would you like a product which:

- **Does not require any special training?**
- **Can test most surfaces ?**
- **Does not require a microbiology laboratory to test the samples?**

If you have answered "YES" to any of these questions, then the STtuve III is the product for you! Available in 4 versions, for the detection of:

- **Salmonella**
- **Vibrio parahaemolyticus**
- **Coliforms/Escherichia coli**
- **Staphylococcus aureus**